

Bacteria Datasheet



City					
Sample Collectors		Date	Start Time	End Time	
Lab Processors				Reminders	
				Rinse bottle 3x and collect 4th	
				Fill to round shoulder of bottle	
				Deliver to lab in < 6 hrs	
				E. coli threshold > 320 cfu/100 ml	
Colony Counter				Incubator IN	Incubator OUT

Site ID	Bottle ID	Sample Vol ml	E. coli Count	Multiply	E. coli cfu/100ml	Comments
-dup		3		33.3		field duplicate
		3		33.3		
		3		33.3		
		3		33.3		
		3		33.3		
		3		33.3		
Distilled water						negative control

Notes: including observations of odor, color, or other visual indicators of pollution